






















# Contents

- **Golden vac** **01**
- **Specimen containers** **14**
- **Petri dish & loops** **26**
- **Nursing plastic ware** **30**
- **Centrifuge tubes & racks** **37**
- **Pipette tips & boxes** **41**
- **Culture plate & pipettes** **45**
- **Test tubes & racks** **51**
- **CUVETTES** **58**
- **CUVETTES** **65**
- **Pharmaceutical packaging material** **68**

GOLDEN VAC

| Category                    | Item                      | Additive                              | Cap Color | Tube material  | Tube cap  | Tube Size (mm)            | Test Item   |
|-----------------------------|---------------------------|---------------------------------------|-----------|----------------|---|---------------------------|---|
| Serum Blood Collection Tube | Plain Tube                | Plain                                 | Red       | Glass/ Plastic |    | 13*75<br>13*100<br>16*100 | Clinical Biochemistry, Immunology and Serology Test |
|                             |                           |                                       |           |                |    | 16*100                    |   |
|                             | Pro-coagulation Tube      | Clot Activator                        | Red       | Glass/ Plastic |    | 13*75<br>13*100<br>16*100 |   |
|                             |                           |                                       |           |                |    | 16*100                    |   |
|                             | Gel & Clot Activator Tube | Gel & Coagulant                       | Yellow    | Glass/ Plastic |    | 13*75<br>13*100<br>16*100 |   |
|                             |                           |                                       |           |                |    | 16*100                    |   |
| Whole Blood Collection Tube | EDTA Tube                 | Sprayed K2EDTA/ Sprayed K3EDTA        | Purple    | Glass/ Plastic |   | 13*75<br>13*100<br>16*100 | Haematology Test (Blood Routine Examination)        |
|                             |                           |                                       |           |                |  | 16*100                    |   |
|                             | ESR Tube                  | 3.8%Sodium Citrate Buffer(0.129mol/L) | Black     | Glass/ Plastic |  | 13*75                     | Erythrocyte Sedimentation Rate Test                 |
|                             |                           |                                       |           | Glass          |  | 8*120                     |   |

GOLDEN VAC

| Category                     | Item               | Additive  | Cap Color | Tube material  | Tube cap  | Tube Size (mm)            | Test Item  |
|------------------------------|--------------------|---|-----------|----------------|---|---------------------------|--|
| Plasma Blood Collection Tube | Coagulation Tube   | 3.2%Sodium Citrate Buffer(0.109mol/L)   | Blue      | Glass/ Plastic |    | 13*75<br>13*100           | Coagulation Function Test                                      |
|                              |                    |   |           |                |    | 10.25*64                  |  |
|                              | Heparin Tube       | Sodium Heparin/ Lithium Heparin   | Green     | Glass/ Plastic |    | 13*75<br>13*100<br>16*100 | Clinical Chemistry for Emergency Treatment,Blood Rheology Test |
|                              |                    |   |           |                |    | 16*100                    |  |
|                              | Gel & Heparin Tube | Gel & Sodium Heparin/ Gel & Lithium Heparin   | Green     | Glass/ Plastic |    | 13*75<br>13*100<br>16*100 | Glucose and Lactate Test                                       |
|                              | Glucose Tube       | Sodium Fluoride & Sodium Heparin/ Sodium Fluoride & EDTA/ Sodium Fluoride & Potassium Oxalate | Grey      | Glass/ Plastic |    | 13*75<br>13*100           |  |
|                              | EDTA & Gel Tube    | Gel & Sprayed K2EDTA/ Gel & Sprayed K3EDTA  | Purple    | Glass/ Plastic |   | 13*75<br>13*100<br>16*100 | Molecular Biology Test (Such as PCR)                           |
|                              |                    |   |           |                |  | 16*100                    |  |


Serum Blood Collection Tube

Golden Vac™ serum blood collection tube promises to provide non-polluting and non-differential original blood specimen for medical test. After centrifugation, serum can be separated from blood cells and fibrin effectively. Serum tube includes 3 kinds of tubes: Plain Tube with red cap, Pro-coagulation Tube with red cap and Gel & Clot Activator Tube with yellow cap.



Plain Tube aims to obtain serum specimen for biochemistry, immunology and serology tests in clinical inspection. Its main merit is that serum separate out of blood cells and fibrin naturally, so its physico-chemical properties are keeping nonpolluted. The intact serum is suitable for a wide range of laboratory operation, routine clinical chemistry, immunology and serology.


| Item   | Specification         |                | Qty/Cs |
|--------|-----------------------|----------------|--------|
| 020A   | Red cap, Plain        | 13*75mm, 2ml   | 1200   |
| 030A   | Red cap, Plain        | 13*75mm, 3ml   | 1200   |
| 040A   | Red cap, Plain        | 13*75mm, 4ml   | 1200   |
| 050A   | Red cap, Plain        | 13*75mm, 5ml   | 1200   |
| 060A   | Red cap, Plain        | 13*100mm, 6ml  | 1200   |
| 070A   | Red cap, Plain        | 13*100mm, 7ml  | 1200   |
| 090A   | Red cap, Plain        | 16*100mm, 9ml  | 1200   |
| 0100A  | Red cap, Plain        | 16*100mm, 10ml | 1200   |
| 090AR  | Rubber stopper, Plain | 16*100mm, 9ml  | 1200   |
| 0100AR | Rubber stopper, Plain | 16*100mm, 10ml | 1200   |



Pro-coagulation Tube aims to obtain serum specimen, which is used for biochemistry and immunology tests. The interior wall is finely coated with coagulant which may promote the blood coagulation after collection. The appropriate amount of coagulant make blood coagulation in a relatively short time, avoiding potential hemolysis problem which caused by too rapid clotting. Finally, transparent serum can be separated after centrifugation.

| Item    | Specification                  |                | Qty/Cs |
|---------|--------------------------------|----------------|--------|
| 020CA   | Red cap, Clot Activator        | 13*75mm, 2ml   | 1200   |
| 030CA   | Red cap, Clot Activator        | 13*75mm, 3ml   | 1200   |
| 040CA   | Red cap, Clot Activator        | 13*75mm, 4ml   | 1200   |
| 050CA   | Red cap, Clot Activator        | 13*75mm, 5ml   | 1200   |
| 060CA   | Red cap, Clot Activator        | 13*100mm, 6ml  | 1200   |
| 070CA   | Red cap, Clot Activator        | 13*100mm, 7ml  | 1200   |
| 090CA   | Red cap, Clot Activator        | 16*100mm, 9ml  | 1200   |
| 0100CA  | Red cap, Clot Activator        | 16*100mm, 10ml | 1200   |
| 090CAR  | Rubber stopper, Clot Activator | 16*100mm, 9ml  | 1200   |
| 0100CAR | Rubber stopper, Clot Activator | 16*100mm, 10ml | 1200   |

Gel & Clot Activator Tube aims to obtain high quality serum specimen for biochemistry, immunology and serology tests in clinical inspection, Inside the tube there is some gel present at the bottom, which is a pure substance and is very stable in physico-chemistry properties. Otherwise, the gel can avoid high-temperature decomposition and small molecules precipitation, hence it greatly reduces the risk of testing equipment clogging.



| Item    | Specification            |                                    | Qty/Cs |
|---------|--------------------------|------------------------------------|--------|
| 040SGC  | Yellow cap, Separate     | Gel+Clot Activator 13*75mm, 4ml    | 1200   |
| 040SGC2 | Yellow cap, Separate     | Gel+Clot Activator 13*100mm, 4ml   | 1200   |
| 050SGC  | Yellow cap, Separate     | Gel+Clot Activator 13*100mm, 5ml   | 1200   |
| 060SGC  | Yellow cap, Separate     | Gel+Clot Activator 13*100mm, 6ml   | 1200   |
| 080SGC  | Yellow cap, Sepearate    | Gel+Clot Activator 16*100mm, 8ml   | 1200   |
| 085SGC  | Yellow cap, Sepearate    | Gel+Clot Activator 16*100mm, 8.5ml | 1200   |
| 080SGCR | Rubber stopper, Separate | Gel+Clot Activator 16*100mm, 8ml   | 1200   |
| 085SGCR | Rubber stopper, Separate | Gel+Clot Activator 16*100mm, 8.5ml | 1200   |

Whole Blood Collection Tube

Golden Vac™ whole blood collection tube promises to provide non-polluting and non-differential original blood specimen for medical test, doesn't require separation of plasma or cells. It includes 2 kinds of tubes: EDTA Tube with Purple cap, ESR Tube with black cap.



EDTA Tube is widely used for clinical haematology test, especially for routine blood examination. Our automatic sampling instrument finely and accurately spray additives onto the interior tube wall, to ensure the mixing speed and the anticoagulation effect. EDTA Tube is very convenient to use for it can be directly put on the analyzer to detect without opening the cap.



| Item     | Specification            |                |  | Qty/Cs |
|----------|--------------------------|----------------|--|--------|
| 020EK3   | Purple cap, EDTA K3,     | 13*75mm, 2ml   |  | 1200   |
| 030EK3   | Purple cap, EDTA K3,     | 13*75mm, 3ml   |  | 1200   |
| 040EK3   | Purple cap, EDTA K3,     | 13*75mm, 4ml   |  | 1200   |
| 050EK3   | Purple cap, EDTA K3,     | 13*75mm, 5ml   |  | 1200   |
| 060EK3   | Purple cap, EDTA K3,     | 13*100mm, 6ml  |  | 1200   |
| 070EK3   | Purple cap, EDTA K3,     | 13*100mm, 7ml  |  | 1200   |
| 080EK3   | Purple cap, EDTA K3,     | 16*100mm, 8ml  |  | 1200   |
| 090EK3   | Purple cap, EDTA K3,     | 16*100mm, 9ml  |  | 1200   |
| 0100EK3  | Purple cap, EDTA K3,     | 16*100mm, 10ml |  | 1200   |
| 080EK3R  | Rubber stopper, EDTA K3, | 16*100mm, 8ml  |  | 1200   |
| 090EK3R  | Rubber stopper, EDTA K3, | 16*100mm, 9ml  |  | 1200   |
| 0100EK3R | Rubber stopper, EDTA K3, | 16*100mm, 10ml |  | 1200   |
| 020EK2   | Purple cap, EDTA K2,     | 13*75mm, 2ml   |  | 1200   |
| 030EK2   | Purple cap, EDTA K2,     | 13*75mm, 3ml   |  | 1200   |
| 040EK2   | Purple cap, EDTA K2,     | 13*75mm, 4ml   |  | 1200   |
| 050EK2   | Purple cap, EDTA K2,     | 13*75mm, 5ml   |  | 1200   |
| 060EK2   | Purple cap, EDTA K2,     | 13*100mm, 6ml  |  | 1200   |
| 070EK2   | Purple cap, EDTA K2,     | 13*100mm, 7ml  |  | 1200   |
| 080EK2   | Purple cap, EDTA K2,     | 16*100mm, 8ml  |  | 1200   |
| 090EK2   | Purple cap, EDTA K2,     | 16*100mm, 9ml  |  | 1200   |
| 0100EK2  | Purple cap, EDTA K2,     | 16*100mm, 10ml |  | 1200   |
| 080EK2R  | Rubber stopper, EDTA K2  | 16*100mm, 8ml  |  | 1200   |
| 090EK2R  | Rubber stopper, EDTA K2  | 16*100mm, 9ml  |  | 1200   |
| 0100EK2R | Rubber stopper, EDTA K2, | 16*100mm, 10ml |  | 1200   |

| Item     | Specification                         |                 |  | Qty/Cs |
|----------|---------------------------------------|-----------------|--|--------|
| 035EK3G  | Purple cap, EDTA K3+Separate Gel,     | 13*75mm, 3.5ml  |  | 1200   |
| 060EK3G  | Purple cap, EDTA K3+Separate Gel,     | 13*100mm, 6ml   |  | 1200   |
| 085EK3RG | Rubber stopper, EDTA K3+Separate Gel, | 16*100mm, 8.5ml |  | 1200   |
| 085EK3G  | Purple cap, EDTA K3+Separate Gel,     | 16*100mm, 8.5ml |  | 1200   |

ESR Tube



ESR Tube is mainly used for Erythrocyte Sedimentation Rate test, it contains buffered Sodium Citrate solution with the concentration of 3.8%(0.129mol/L).The mixing ratio of blood to Sodium Citrate is 4:1.It is suitable for the ESR tests for subjects of different gender and ages, with highly related coefficient comparing with the Westertren method..

| Item  | Specification  |                |  | Qty/Cs |
|-------|----------------|----------------|--|--------|
| 016SC | Black cap, ESR | 13*75mm, 1.6ml |  | 1200   |
| 024SC | Black cap, ESR | 13*75mm, 2.4ml |  | 1200   |
| 032SC | Black cap, ESR | 13*75mm, 3.2ml |  | 1200   |

ESR Tube



ESR Detector

Condition: Within 6 hours to finish the testing.  
Lab temperature: 8 to 25℃

Instructions:  
1. Put the ESR detector on the level place.  
2. Using the ESR tube to collect the blood and mixing well.  
3. Put the tube upright on the ESR detector.The tube blue graduation line flush with rack top.  
4. Wait 60mins, then read the data that the level of blood go down.  
Advantage:


ESR Fast Dectector have a Particular organic glass structure that make the ESR tubes fixed fast and simply, promise it the perpendicular degree of the level direction, and can conveniently adj ust the horizon of ESR tubes samples.

| Item    | Specification             |                     |  | Qty/Cs |
|---------|---------------------------|---------------------|--|--------|
| 016ESR  | Rubber stopper, ESR Tube  | 8*120mm, 1.6ml      |  | 1200   |
| 0128ESR | Rubber stopper, ESR Tube  | 8*120mm, 1.28ml     |  | 1200   |
| 020A    | ESR Fast Detector, for    | 8X120mm,1.6ml tube  |  | 10     |
| 020B    | ESR Fast Detector, for    | 8X120mm,1.28ml tube |  | 10     |
| 016ESRN | Rubber Stopper, ESR Tube, | 9x120mm, 1.6ml      |  | 1200   |



Plasma Blood Collection Tube

Golden Vac™ plasma blood collection tube promises to provide non-polluting and non-differential original blood specimen for medical test. After centrifugation, plasma can be separated from blood cells effectively. Plasma tube includes 5 kinds of tubes: Clotting Tube with blue cap, Heparin Tube and Gel & Heparin Tube with green cap, Compound of Sodium Fluoride and Anticoagulant Tube with grey cap, Gel & EDTA Tube with lavender cap.





Coagulation Tube aims to obtain plasma specimen for coagulation function test, it contains buffered Sodium Citrate solution with the concentration of 3.2%(0.109mol/L).The mixing ratio of blood to Sodium Citrate is 9:1. This anticoagulant is sustained and stable, it strictly control the microbial limit, Making sure the APTT test under the excellent Condition.

| Item   | Specification   | Qty/Cs |
|--------|---|--------|
| 018SC  | Blue cap, Coagulation, 3.2% Sodium citrate 13*75mm, 1.8ml | 1200   |
| 027SC  | Blue cap, Coagulation, 3.2% Sodium citrate 13*75mm, 2.7ml | 1200   |
| 036SC  | Blue cap, Coagulation, 3.2% Sodium citrate 13*75mm, 3.6ml | 1200   |
| 045SC  | Blue cap, Coagulation, 3.2% Sodium citrate 13*75mm, 4.5ml | 1200   |
| 018SCD | Blue cap, Coagulation, 3.2% Double wall 13*75mm, 1.8ml    | 1200   |
| 027SCD | Blue cap, Coagulation, 3.2% Double wall 13*75mm, 2.7ml    | 1200   |



| Item    | Specification  | Qty/Cs |
|---------|--|--------|
| 018SC-F | Blue cap, Glass tube,Coagulation, 10.25*64mm(flat bottom), 1.8ml | 1000   |
| 027SC-F | Blue cap, Glass tube,Coagulation, 10.25*64mm(flat bottom), 2.7ml | 1000   |



(Aluminium Plastic bag pack for each rack)



Heparin Tube aims to obtain plasma specimen, which is used not only for emergency biochemistry tests but also for blood rheology test. Inside the tube wall, where is well sprayed with lithium heparin or sodium heparin which is a kind of anticoagulant activating antithrombins to block the coagulation cascade reaction.

| Item    | Specification                                   | Qty/Cs |
|---------|---|--------|
| 030SH   | Green cap, Heparin Sodium, 13*75mm, 3ml         | 1200   |
| 040SH   | Green cap, Heparin Sodium, 13*75mm, 4ml         | 1200   |
| 050SH   | Green cap, Heparin Sodium, 13*75mm, 5ml         | 1200   |
| 060SH   | Green cap, Heparin Sodium, 13*100mm, 6ml        | 1200   |
| 070SH   | Green cap, Heparin Sodium, 13*100mm, 7ml        | 1200   |
| 0100SH  | Green cap, Heparin Sodium, 16*100mm, 10ml       | 1200   |
| 030LH   | Green cap, Heparin Lithium, 13*75mm, 3ml        | 1200   |
| 040LH   | Green cap, Heparin Lithium, 13*75mm, 4ml        | 1200   |
| 050LH   | Green cap, Heparin Lithium, 13*75mm, 5ml        | 1200   |
| 060LH   | Green cap, Heparin Lithium, 13*100mm, 6ml       | 1200   |
| 070LH   | Green cap, Heparin Lithium, 13*100mm, 7ml       | 1200   |
| 0100LH  | Green cap, Heparin Lithium, 16*100mm, 10ml      | 1200   |
| 0100LHR | Rubber stopper, Heparin Lithium, 16*100mm, 10ml | 1200   |



Gel & Heparin Tube is a kind of anticoagulation tube, aims to obtain high quality plasma specimen for fast emergency biochemistry test. Lying at the bottom of the tube, there is some gel for separating plasma and blood cells, preventing substance from exchanging and keeping the obtained plasma from being blood cell polluted, thus the accuracy of the results will be improved. Otherwise, the specimen could be put into the equipment directly without transferring and kept stable in cold storage for 4hs, so that well meets the quality requirement of retest.

| Item   | Specification                                   | Qty/Cs |
|--------|---|--------|
| 030SHG | Green cap, Heparin Sodium+Gel, 13*75mm, 3ml     | 1200   |
| 040SHG | Green cap, Heparin Sodium+Gel, 13*75mm, 4ml     | 1200   |
| 050SHG | Green cap, Heparin Sodium+Gel, 13*100mm, 5ml    | 1200   |
| 085SHG | Green cap, Heparin Sodium+Gel, 16*100mm, 8.5ml  | 1200   |
| 030LHG | Green cap, Heparin Lithium+Gel, 13*75mm, 3ml    | 1200   |
| 040LHG | Green cap, Heparin Lithium+Gel, 13*75mm, 4ml    | 1200   |
| 050LHG | Green cap, Heparin Lithium+Gel, 13*100mm, 5ml   | 1200   |
| 085LHG | Green cap, Heparin Lithium+Gel, 16*100mm, 8.5ml | 1200   |

## GOLDEN VAC

## Glucose Tube



Compound of Sodium Fluoride and Anticoagulant Tube aims to obtain highly stable plasma for blood glucose testing .Sodium Fluoride is a kind of inhibitor of glucose degradation ,it has the function of repressing glucose metabolism and degradation .Optional Anticoagulant are Potassium Oxalate ,EDTA and Heparin ,they could prevent blood coagulation and avoid blood hemolysis .This particular compound additive realise that blood glucose specimen could be cold storage for a long time.

| Item   | Specification   | Qty/Cs |
|--------|---|--------|
| 020SFH | Gray cap, Sodium Fluoride+Heparin Sodium, 13*75mm, 2ml  | 1200   |
| 030SFH | Gray cap, Sodium Fluoride+Heparin Sodium, 13*75mm, 3ml  | 1200   |
| 040SFH | Gray cap, Sodium Fluoride+Heparin Sodium, 13*75mm, 4ml  | 1200   |
| 050SFH | Gray cap, Sodium Fluoride+Heparin sodium, 13*75mm, 5ml  | 1200   |
| 020SNA | Gray cap, Sodium Fluoride+EDTA Naz, 13*75mm, 2ml  | 1200   |
| 030SNA | Gray cap, Sodium Fluoride+EDTA Naz, 13*75mm, 3ml  | 1200   |
| 050SNA | Gray cap, Sodium Fluoride+EDTA Naz, 13*75mm, 5ml  | 1200   |
| 020SP  | Gray cap, Sodium Fluoride+Potassium Oxalate, 13*75mm, 2ml   | 1200   |
| 030SP  | Gray cap, Sodium Fluoride+Potassium Oxalate, 13*75mm, 3ml   | 1200   |
| 040SP  | Gray cap, Sodium Fluoride+Potassium Oxalate, 13*75mm, 4ml   | 1200   |
| 050SP  | Gray cap, Sodium Fluoride+Potassium Oxalate, 13*75mm, 5ml   | 1200   |
| 6001   | Disposable ESR Pipette with prefilled (3.2% Trisodium Citrate) and capped filling tube 100pcs/box*10boxes |        |

## Non Vacuum Blood Tube



| Item    | Specification   | Qty/Cs |
|---------|---|--------|
| N050EK3 | 5ml EDTA K3 Non-vacuum blood collection tube, Green cap             | 1000   |
| N050EK2 | 5ml EDTA K2 Non-vacuum blood collection tube, Green cap             | 1000   |
| N050A   | 5ml Plain Non-vacuum blood collection tube, Red cap                 | 1000   |
| N050CA  | 5ml Clot Activator Non-vacuum blood collection tube, Red cap        | 1000   |
| N050SP  | 5ml Fluoride Oxalate Non-vacuum blood collection tube, Yellow cap   | 1000   |
| N050LH  | 5ml Lithium Heparin Non-vacuum blood collection tube, Dark Blue cap | 1000   |



## Micro Collection Tube



Micro Collection tube is mainly used for peripheral blood collection. It is specially suitable for small children and geriatric patient. The inner surface has an excellent hemorepellency where it is necessary for a good preparation of high quality blood sample.

## Capillary Blood Collection Tube



## GOLDEN VAC

| Item     | Specification  | Qty/Cs |
|----------|--|--------|
| 005A     | Micro Collection tube, Plain, 0.5ml                                    | 1000   |
| 005CA    | Micro Collection tube, Clot Activator 0.5ml                            | 1000   |
| 005EK3   | Micro Collection tube, EDTA K3, 0.5ml                                  | 1000   |
| 005EK2   | Micro Collection tube, EDTA K2, 0.5ml                                  | 1000   |
| 005LH    | Micro Collection tube, Heparin Lithium, 0.5ml                          | 1000   |
| 005SH    | Micro Collection tube, Heparin Sodium, 0.5ml                           | 1000   |
| 005SFH   | Micro Collection tube, Heparin Sodium +Sodim Fluoride, 0.5ml           | 1000   |
| 005SGC   | Micro Collection tube, Gel+Clot Activator, 0.5ml                       | 1000   |
| NM005A   | Micro tube (High Type) 0.5ml Plain 0.5ml                               | 1000   |
| NM005CA  | Micro tube (High Type) 0.5ml Clot activator 0.5ml                      | 1000   |
| NM005EK3 | Micro tube (High Type) 0.5ml EDTAK3 0.5ml                              | 1000   |
| NM005EK2 | Micro tube (High Type) 0.5ml EDTAK2 0.5ml                              | 1000   |
| NM005LH  | Micro tube (High Type) 0.5ml, Gel 0.5ml                                | 1000   |
| NM005SGC | Micro tube (High Type) 0.5ml, Gel 0.5ml                                | 1000   |
| NM005SFH | Micro tube (High Type) Glucose(Heparin sodium,/ Sodium Fluoride) 0.5ml | 1000   |




| Item    | Specification   | Qty/Cs |
|---------|---|--------|
| C005A   | Capillary Blood Collection tube, Plain, 0.5ml                         | 1000   |
| C005CA  | Capillary Blood Collection tube, Clot Activator, 0.5ml                | 1000   |
| C005EK3 | Capillary Blood Collection tube, EDTA K3, 0.5ml                       | 1000   |
| C005EK2 | Capillary Blood Collection tube, EDTA K2, 0.5ml                       | 1000   |
| C005LH  | Capillary Blood Collection tube, Heparin Lithium, 0.5ml               | 1000   |
| C005SH  | Capillary Blood Collection tube, Heparin Sodium, 0.5ml                | 1000   |
| C005SFH | Capillary Blood Collection tube, Heparin Sodium+Sodim Fluoride, 0.5ml | 1000   |
| C005SGC | Capillary Blood Collection tube, Gel+Clot Activator, 0.5ml            | 1000   |



SHARPS CONTAINER

Sharps Container

FDA 510(K)



For sharp container of 3L, 5L, 7L, two kinds of caps are both available.

| Item      | Specification                           | Qty/Cs |
|-----------|---|--------|
| GDC-01-F  | Sharp container, Flip+Ype closure 1L    | 60     |
| GDC-03-F  | Sharp container, Flip+Ype closure 3L    | 48     |
| GDC-03-X  | Sharp container, Rotary+Ype closure 3L  | 48     |
| GDC-05-F  | Sharp container, Flip+Ype closure 5L    | 32     |
| GDC-05-X  | Sharp container, Rotary+Ype closure 5L  | 32     |
| GDC-07-F  | Sharp container, Flip+Ype closure 7L    | 24     |
| GDC-07-X  | Sharp container, Rotary+Ype closure 7L  | 24     |
| GDC-013-X | Sharp container, Rotary+Ype closure 13L | 20     |

SHARPS CONTAINER



| Item   | Description of Goods                     | Cs size       | Qty/Cs     |
|--------|--|---------------|------------|
| SC-01Q | Sharps Container 1 Quart                 | 108*88*170mm  | 60pcs/ctn  |
| SC-02Q | Sharps Container 2 Quart                 | 108*88*220mm  | 60pcs/ctn  |
| SC-05Q | Sharps Container 5 Quarts                | 273*115*282mm | 20pcs/ctn  |
| SC-01G | Sharps Container 1 Gallons               | 258*175*142mm | 24pcs/ctn  |
| SC-02G | Sharps Container 2 Gallons               | 274*175*254mm | 24pcs/ctn  |
| SC-06G | Sharps Container 6 Gallons               | 378*239*391mm | 12pcs/ctn  |
| SCT50  | Sharps shafe 6.4X1.22 in(16.17X3.1cm)Lxw | 500           | 500pcs/ctn |




GOLDEN VAC



Needle & Holder

| Item   | Specification             | Qty/Cs |
|--------|---------------------------|--------|
| N18G   | Needle 18Gx1 1/2"         | 4000   |
| N20G   | Needle 20Gx1 1/2"         | 4000   |
| N21G   | Needle 21Gx1 1/2"         | 4000   |
| N22G   | Needle 22Gx1 1/4"         | 4000   |
| BF22G  | Butterfly Needle 22Gx3/4  | 2000   |
| N21GFB | Flash Back Needle 21Gx1"  | 5000   |
| N22GFB | Flash Back Needle 22Gx1"  | 5000   |
| H      | Holder                    | 2000   |
| H1     | Holder With Safty Cover   | 2000   |
| H2     | Holder With Quick Release | 2000   |

| Item | Description of Goods   | Specification                     | Packing       |
|------|------------------------|-----------------------------------|---------------|
| B-01 | blood collect tube box | 13*75mm                           | 50pcs/in bulk |
| B-02 | blood collect tube box | 13*100mm                          | 50pcs/in bulk |
| S-01 | specimen storage box   | 13*100mm ( pillar base )100 holes | 28pcs/cs      |
| S-02 | specimen storage box   | 13*75mm ( pillar base )100 holes  |               |







Specimen containers








SPECIMEN CONTAINERS



SPECIMEN CONTAINERS




Stool specimen collection

It is for the collection of stool specimen and is clean, convenient and safe. It can form the closed collection and transportation of stool specimen, which effectively prevents the drying, spill and leakage of stool specimen and is good for specimen marking, bulk delivering, and Sampling during specimen testing.

| Item    | Description of Goods                     | Quantity (pcs) | Packing       |
|---------|--|----------------|---------------|
| K1014   | 5ml stool collection tube(with spoon)    | 2000           | 25pcs*80bags  |
| K1014-1 | 5ml stool collection tube(without spoon) | 2000           | 25pcs*80bags  |
| K1018-2 | Stool container with spoon 20ml          | 650            | 25pcs*26bags  |
| K1018-3 | Stool container with spoon 20ml          | 650            | in bulk       |
| K1011   | 20ml stool container                     | 1000           | 100pcs*10bags |
| K1012   | 25ml culture bottle                      | 650            | 25pcs*26bags  |
| K1012-1 | 25ml stool container                     | 650            | 25pcs*26bags  |

SPECIMEN CONTAINERS



| Item    | Description of Goods                   | Quantity <sub>(pcs)</sub> | Packing       |
|---------|--|---------------------------|---------------|
| K1017-2 | 30ml universal container with spoon PP | 650                       | 25pcs*26bags  |
| K1017-3 | 30ml universal container with spoon PS | 650                       | 25pcs*26bags  |
| K1013   | 60ml stool container                   | 500                       | in bulk       |
| K1013-1 | 60ml stool container                   | 500                       | in bulk       |
| K1013-2 | Stool stick                            | 5000                      | 100pcs*50bags |

SPECIMEN CONTAINERS



Sputum specimen collection

It is for the collection of sputum specimen and can form the closed collection and transportation of sputum specimen. It is clean, convenient and safe.


| Item    | Description of Goods                   | Quantity <sub>(pcs)</sub> | Packing      |
|---------|--|---------------------------|--------------|
| K1015   | 30ml sputum container                  | 500                       | 1pc*500bags  |
| K1018-1 | 20ml sputum container(NEW)             | 1000                      | in bulk      |
| K1017   | 30ml universal container               | 650                       |              |
| K1017-1 | 30ml sputum container (skirted bottom) | 650                       | 25pcs*26bags |

SPECIMEN CONTAINERS



| Item  | Description of Goods  | Quantity <sub>(pcs)</sub> | Packing        |
|-------|-----------------------|---------------------------|----------------|
| K1018 | 30ml sputum container | 1000                      | 10pcs*100bags  |
| K1016 | 40ml sputum container | 500                       | 5pcs*100bags   |
| K1019 | Sputum cassette       | 2000                      | 100pcs*10bags  |
| K1023 | Mucus trap,40cc       | 50                        | 1pc/bag*50bags |
| K1024 | Mucus trap,80cc       | 50                        | 1pc/bag*50bags |

SPECIMEN CONTAINERS



Urine specimen collection

It is for the collection of urine specimen. It is made of medical plastic and designed in humanization with thick and even wall. It can form the closed collection, transportation and testing of urine, which has solved a series of problems, such as urine specimen submission, supporting the use of computers, pollution and so on.


- The product use is clean, convenient and safe, which provides a guarantee for the normalization and standardization of urine analysis from the source.

| Item  | Description of Goods                        | Quantity <sub>(pcs)</sub> | Packing    |
|-------|---|---------------------------|------------|
| K1009 | Soft urine cup white                        | 10000                     | in bulk    |
| K1010 | Soft urine cup transparent                  | 10000                     | in bulk    |
| K1025 | 30ml urine container with mouth on cap,30ml | 500                       | 1pc*500pcs |

SPECIMEN CONTAINERS




SPECIMEN CONTAINERS



| Item    | Description of Goods           | Quantity (pcs) | Packing                                       |
|---------|--------------------------------|----------------|---|
| H1016   | 60ml specimen container        | 500            | 1pc*500pcs                                    |
| H1016-1 | 60ml specimen container, PS    | 500            | 1pc*500pcs                                    |
| H1016-2 | 90ml specimen container        | 500            | 5pcs*100bags                                  |
| H1016-3 | 40ml specimen container        | 500            | 1pc*500pcs                                    |
| H1016-5 | 40ml renovate urine container  | 500            | 50pcs/ctn                                     |
| H1015   | 100ml-120ml specimen container | 500/1000       | Individual Packing/ 20pcs*50bags Bulk Packing |
| H1021   | 120ml-140ml Urine container    | 500            | 1pc*500pcs                                    |



SPECIMEN CONTAINERS



| Item     | Description of Goods  | Quantity <sub>(pcs)</sub> | Packing       |
|----------|---|---------------------------|---------------|
| K1026    | Urine Transfer Device for Urine Tube, 100ml                 | 250                       | 1pc*250pcs    |
| GD0105UA | Vacuum Urine Tube, Yellow cap No Additive, 16*100mm, 10.5ml | 1200                      | 1pc*1200pcs   |
| GD0105UB | Vacuum Urine Tube, Yellow cap Boric Acid, 16*100mm, 10.5ml  | 1200                      | 1pc*1200pcs   |
| H1026    | Urine Storage Container for 24hours, 3000ml                 | 50                        |               |
| H1027    | Urine tube, 16X85mm, 10ml                                   | 2000                      | 200pcs*10bags |

TRANSPORT SWABS WITH MEDIUM



Transport Swabs, Amies, Snappable Plastic Stick + Viscose head, Swab length:150mm. Tubes in PP material, Dia.13\*150mm, labelled, Gamma Sterile. Individual Packing. Blue Cap



Transport Swabs, Amies with Charcoal, Snappable Plastic Stick + Viscose head, Swab length:150mm. Tubes in PP material, Dia.13\*150mm, labelled, Gamma Sterile. Individual Packing. White Cap



Transport Swabs, Cary Blair, Snappable Plastic Stick+Viscose head, Swab length: 150mm. Tubes in PP material, Dia.13\*150mm, labelled, Gamma Sterile. Individual Packing. Red Cap



Transport Swabs, Snappable Plastic Stick+ Viscose head, Swab length:150mm. Tubes in PP material, Dia.13\*150mm, labelled, Gamma Sterile. Individual Packing. Black Cap



Transport Swabs, Stuart with Charcoal, Snappable Plastic Stick+Viscose head, Swab length:150mm. Tubes in PP material, Dia.13\*150mm, labelled, Gamma Sterile. Individual Packing. Black Cap

| Item   | Description of Goods                              | Quantity <sub>(pcs)</sub> | Packing       |
|--------|---|---------------------------|---------------|
| G10261 | Transport swab with medium (Amies)                | 600                       | 50pcs*12boxes |
| G10262 | Transport swab with medium (Amies with charcoal)  | 600                       | 50pcs*12boxes |
| G10263 | Transport swab with medium (Cary Blair)           | 600                       | 50pcs*12boxes |
| G10264 | Transport swab with medium (Stuart)               | 600                       | 50pcs*12boxes |
| G10265 | Transport swab with medium (Stuart with charcoal) | 600                       | 50pcs*12boxes |

TRANSPORT SWABS WITH MEDIUM



| Item  | Description of Goods            | Quantity (pcs) | Packing       |
|-------|---------------------------------|----------------|---------------|
| G1022 | 13*150mm, tube with Swab        | 1800           | 150pcs*12bags |
| G1023 | 12*75mm, tube with Swab         | 2000           |               |
| G1028 | disposable transport sample box |                |               |



Petri dish& loops

SYSTEM

Raycaremed Petri dish series are made of imported polystyrene with high quality as raw material. Its smoothness and transparent finish fully achieve the International standard of similar product. It is produced in a 100 thousand stage purification plant where there is no heat source and sterilized by ethylene oxide. Adopt the internationally popular easy tear package, which is simple and convenient and reduce the pollution opportunities in the use of products. Each batch of products is evaluated by growth status observation, test of cell growth rate, test of cell adherence strength, pyrogenicity test and sterility test to ensure product quality In comparison with similar quality products in the world.

petri dish series is suitable for cell, bacteria culture and drug sensitivity test.

PETRI DISH & LOOPS



| Item  | Description of Goods          | Quantity <sub>(pcs)</sub> | Packing            |
|-------|-------------------------------|---------------------------|--------------------|
| K1001 | 15cm Petri dish               | 180                       | Individual Packing |
| K1002 | 9cm Petri dish $\phi$ 9*2.0cm | 500                       | 10pcs*50bags       |
| K1020 | 10ul Sampling Loop            | 3000                      | 10pcs/bag*300bags  |
| K1021 | 1ul Sampling Loop             | 3000                      | 10pcs/bag*300bags  |
| K1022 | Needle Type Sampling Loop     | 3000                      | 10pcs/bag*300bags  |



PETRI DISH




● K1003


Facilitate the differentiation of bacterial culture programs


● K1004

There are three points on the cover to prevent the accumulation of bacterial culture specimens.

| Item            | Description of Goods           | Quantity <sub>(pcs)</sub> | Packing          |
|-----------------|--------------------------------|---------------------------|------------------|
| K1003           | 90X15mm Two-Compartment        | 500                       | 10pcs/bag*50bags |
| K1004-J         | 90X15mm Petri Dish             | 500                       | 10pcs/bag*50bags |
| K1004-J-3 Vents | 90X15mm 3Vents Petri Dish      | 500                       | 10pcs/bag*50bags |
| K1004-M         | 90X15mm for machine Petri Dish | 500                       | 10pcs/bag*50bags |
| K1004-T         | 90X15mm 3Vents Petri Dish      | 500                       | 10pcs/bag*50bags |





 **K1006** The split format Petri dish is used for bacterial cultures that require the separation of curved surface. The mesh distance is 10mmX10mm. The numbers and alphabetic numbers are convenient for location, identification and counting.

| Item    | Description of Goods                                | Quantity (pcs) | Packing           |
|---------|---|----------------|-------------------|
| K1005   | 70x15mm Petri dish                                  | 600            | 10pcs/bag*60bags  |
| K1006   | 65x15mm Petri dish                                  | 650            | 10pcs/bag*65bags  |
| K1007   | 55x15mm Petri dish                                  | 1000           | 10pcs/bag*100bags |
| K1008   | 30x15 Petri dish                                    | 2000           | 10pcs/bag*200bags |
| K1009   | 14cm Petri Dishes                                   | 180            | 10pcs*18bags      |
| K1023   | Hook match with mucus trap<br>inner-Φ6.5mm out-Φ9mm | 1000           |                   |
| K1032-L | 29cm Cell Scraper                                   | 100            | Individual        |
| K1032-S | 25cm Cell Scraper                                   | 100            | Individual        |
| K1033   | 18cm Cell Lifter                                    | 100            | Individual        |




# Nursing plastic ware

SYSTEM

Vaginal Speculum got FDA approved without any defect. The products are sold well in the United States and other countries. The product design is delicate, humanized and simple to operate. A variety of different specifications of Dressing Cups are suitable for all ages, all kinds of medicine, convenient and clean. All kinds of PP plastic tweezers adopt connecting lever design, more excellent performance, can meet different working requirements.



SPECULUM



| Item    | Description of Goods                      | Quantity <sub>(pcs)</sub> | Packing          |
|---------|---|---------------------------|------------------|
| P1001   | Vaginal Speculum Large,Middle,Small       | 100                       | 10pcs/bag*10bags |
| P1002-L | Vaginal speculum,Lighe souce style,large  | 80                        | 8pcs/bag*10bags  |
| P1002-M | Vaginal speculum,Lighe souce style,Medium | 100                       | 10pcs/bag*10bags |
| P1002-S | Vaginal speculum,Lighe souce style,Small  | 100                       | 10pcs/bag*10bags |
| P1030   | Anus speculum, Large,Middle,Small         | 100                       | 25pcs/bag*4bags  |
| P1039   | Anus speculum                             | 100                       | 25pcs/bag*4bags  |



FORCEPS



| Item          | Description of Goods                     | Quantity <sub>(pcs)</sub> | Packing           |
|---------------|--|---------------------------|-------------------|
| P1012 a/b/c/e | Plastic Tweezers (sterile )              | 2500                      | 250pcs/bag*10bags |
| P1012-1       | Fine-tipped Tweezers(sterile ) 10cm long | 2500                      | 250pcs/bag*12bags |
| P1013         | Tweezers with tooth(sterile )            | 2500                      | 250pcs/bag*10bags |
| P1014         | Long Tweezers tooth(sterile )            | 2000                      | 200pcs/bag*10bags |
| P1015         | Tweezers without tooth(sterile )         | 2500                      | 250pcs/bag*10bags |
| P1016         | Blue Tweezers (sterile )                 | 2000                      | 200pcs/bag*10bags |



SPONGE HOLDER



33

SPONGE HOLDER

| Item  | Description of Goods           | Quantity <sub>(pcs)</sub> | Packing           |
|-------|--------------------------------|---------------------------|-------------------|
| P1020 | Sponge Holder 23cm Long        | 500                       | 50pcs/bag*10bags  |
| P1021 | Sponge Holder 18cm Long        | 500                       | 50pcs/bag*10bags  |
| P1022 | Sponge Holder 11cm Long        | 1200                      | 100pcs/bag*12bags |
| P1023 | Sponge Holder 10cm Long        | 1200                      | 100pcs/bag*12bags |
| P1024 | Sponge Holder 14cm Long clamp  | 1200                      | 100pcs/bag*12bags |
| P1025 | Sponge Holder 14cm Short clamp | 1200                      | 100pcs/bag*12bags |




NURSING PLASTIC WARE




34

NURSING PLASTIC WARE


| Item  | Description of Goods | Quantity <sub>(pcs)</sub> | Packing                |
|-------|----------------------|---------------------------|------------------------|
| P1004 | Scissors,Blue        | 1200                      | 100pcs/bag*12bags      |
| P1033 | Scissors,Blue        | 1000                      | 100pcs/bag*10bags      |
| P1038 | Tweezres 13cm long   | 2000                      | 200pcs/bag*15bags      |
| K1034 | 2.5mm Ear Speculum   | 8500                      | 34pcs/bag*25bag*10bags |
| k1035 | 4.25mm Ear Speculum  | 8500                      | 34pcs/bag*25bag*10bags |



NURSING PLASTIC WARE



| Item    | Description of Goods                                | Quantity <sub>(pcs)</sub> | Packing            |
|---------|---|---------------------------|--------------------|
| P1008   | Tongue Depressor                                    | 2800                      | 200pcs/bag*14bags  |
| P1010   | Mouth piece for Endoscope Large                     | 1000                      | 100pcs/bag*10bags  |
| P1011   | Mouth piece for Endoscope Medium                    | 1000                      | 100pcs/bag*10bags  |
| P1017   | Bite-block-Large                                    | 2200                      | 200pcs/bag*11bags  |
| P1017-1 | Bite-block-Small                                    | 2700                      | 300pcs/bag*9bags   |
| P1018-1 | Mouth piece for Endoscope with elastic(Large,Small) | 900                       | Individual Packing |



NURSING PLASTIC WARE



| Item    | Description of Goods               | Quantity <sub>(pcs)</sub> | Packing                       |
|---------|------------------------------------|---------------------------|-------------------------------|
| P1005   | 3*1Three-color medicine cup        | 500                       | 25pcs/bag*20bags              |
| P1007   | 3*1white medicine cup              | 500                       | 25pcs/bag*20bags              |
| P1003   | Transparent dressing bowl(Sterile) | 1500                      | 2pcs/packet*50pcs/bags*15bags |
| P1004-1 | Dressing basin(soft)               | 1000                      | 2pcs/packet*50pcs/bags*10bags |
| P1004-2 | Dressing basin(hard )              | 80                        | 2pcs/bag*40bags               |





# Centrifuge tubes & Racks

## SYSTEM

Centrifuge tube series are made of polypropylene with high quality. The thickness of the tube wall is uniform. The tubes are all marked with scales so that the volume of concentrated liquid can be directly read out, and make it easier to recycle concentrated liquid with straw. The unique sealing cover prevents liquid leakage, and can be sterilized in high temperature and pressure. Match with a variety of centrifuge tube holders.

## CENTRIFUGE TUBES & RACKS




| Item     | Description of Goods                         | Quantity(pcs) | Packing      |
|----------|--|---------------|--------------|
| D1001BB  | 50ml conical with black graduation blue cap  | 500           | 50pcs*10bags |
| D1001BW  | 50ml conical with white graduation blue cap  | 500           | 50pcs*10bags |
| D1001GB  | 50ml conical with black graduation green cap | 500           | 50pcs*10bags |
| D1001GW  | 50ml conical with white graduation green cap | 500           | 50pcs*10bags |
| D1002WW  | 50ml self-standing white cap                 | 500           | 50pcs*10bags |
| D1002BB  | 50ml self-standing blue cap                  | 500           | 50pcs*10bags |
| D1003BB  | 15ml conical with black graduation blue cap  | 1000          | 50pcs*20bags |
| D1003BW  | 15ml conical with white graduation blue cap  | 1000          | 50pcs*20bags |
| D1003GB  | 15ml conical with black graduation green cap | 1000          | 50pcs*20bags |
| D1003GW  | 15ml conical with white graduation green cap | 1000          | 50pcs*20bags |
| D1003-PS | 15ml conical with graduation                 | 1000          | 50pcs*20bags |





CENTRIFUGE TUBES & RACKS



| Item  | Description of Goods    | Quantity <sub>(pcs)</sub> | Packing       |
|-------|-------------------------|---------------------------|---------------|
| D1004 | 50ml                    | 1000                      | 50pcs*20bags  |
| D1005 | 15ml conical            | 1500                      | 150pcs*10bags |
| D1006 | 10ml                    | 2000                      | 200pcs*10bags |
| D1008 | 5ml*cap                 | 3000                      | 300pcs*10bags |
| D1013 | 5ml+cap.                | 2000                      | 200pcs*10bags |
| D1015 | 1.8ml                   | 5000                      | 500pcs*10bags |
| D1017 | 1.5 ml linked green cap | 6000                      | 500pcs*12bags |
| D1018 | 1.5ml linked red cap    | 6000                      | 500pcs*12bags |




⦿ D1017 / D1018

It is made of polypropylene with high quality. The thickness of the tube wall is uniform. The connecting pipe cover is provided with a silica gel sealing ring, which has good sealing performance and can prevent leakage pollution during evaporation and transportation.

The color of the cap can be increased according to the needs of the laboratory project.

The thread cap is easy to screw and open.

CENTRIFUGES TUBE & RACKS



| Item  | Description of Goods  | Quantity <sub>(pcs)</sub> | Packing        |
|-------|---|---------------------------|----------------|
| D1009 | 2ml   | 5000                      | 500pcs*10bags  |
| D1010 | 1.5ml with transparent cap  | 6000                      | 500pcs*12bags  |
| D1011 | 0.5ml   | 12000                     | 1000pcs*12bags |
| D1012 | 0.2ml   | 30000                     | 1000pcs*30bags |
| D1021 | 98wells centrifugation tube rack                                      | 40                        | 1pc*40bags     |
| D1022 | 72wells centrifugation tube rack                                      | 40                        | 1pc*40bags     |
| D1023 | 81wells centrifugation tube rack                                      | 40                        | 1pc*40bags     |
| D1024 | two purpose tube rack   | 50                        | 1pc*40bags     |
| D1027 | two-purpose centrifuge tube rack<br>(front 1.5ml, reverse side 0.5ml) | 50                        | 1pc*50bags     |
| H1013 | centrifugation rack orange, for 50ml , 15ml                           | 50                        | 1pc*50bags     |




# Pipette tips & Boxes


## SYSTEM

The special Raycaremed pipette tip is used to ensure the accuracy and precision of the experiment. The precision pipette tip provides optimum reproducibility and accuracy. Standard tip and tip with filter can be disinfected and sterilized in high pressure. The unique filter element silicified tip does not hang the residual liquid, which guarantees the accuracy of liquid transfer to the utmost. It almost meets and fits the requirements of most commercially available pipettes. It is also equipped with a variety of tip boxes.


## PIPETTE TIPS & BOXES




| Item  | Description of Goods                        | Quantity <sub>(pcs)</sub> | Packing            |
|-------|---|---------------------------|--------------------|
| E1001 | φ 4*31mm                                    | 50000                     | 500pcs/bag*100bags |
| E1002 | φ 6*50mm                                    | 50000                     | 500pcs/bag*100bags |
| E1003 | φ 8*71mm                                    | 10000                     | 100pcs/bag*100bags |
| E1004 | 250ul φ 4*49mm                              | 40000                     | 500pcs/bag*80bags  |
| E1005 | 250ul φ 6*50mm                              | 50000                     | 500pcs/bag*100bags |
| E1006 | 1000ul φ 8*71mm                             | 10000                     | 100pcs/bag*100bags |
| E1007 | 1000ul Fit For Eppendorf                    | 10000                     | 100pcs/bag*100bags |
| E1008 | 300ul Fit For Eppendorf                     | 30000                     | 1000pcs/bag*30bags |
| E1009 | 0.1–10ul φ 4*49mm Fit For Eppendorf         | 60000                     | 1000pcs/bag*60bags |
| E1010 | 0.10–10ul φ 5*32mm Fit For Gilson           | 60000                     | 1000pcs/bag*60bags |
| E1011 | 200ul φ 6*50mm ( Yellow) Fit For Gilson     | 50000                     | 500pcs/bag*100bags |
| E1012 | 200ul φ 6*50mm ( Graduation) Fit For Gilson | 30000                     | 1000pcs/bag*30bags |
| E1013 | 1000ul φ 8*71mm Fit For Gilson              | 10000                     | 100pcs/bag*100bags |




PIPETTE TIPS & BOXES



| Item    | Description of Goods                       | Quantity <sub>(pcs)</sub> | Packing            |
|---------|--|---------------------------|--------------------|
| E1014   | 0.10–10ul ϕ 5*32mmFit For Finland          | 60000                     | 1000pcs/bag*60bags |
| E1015   | 40–200ul0.50–10ul ϕ 6*50mm Fit For Finland | 30000                     | 1000pcs/bag*30bags |
| E1016   | 250ul ϕ 6*52mmFit For Finland              | 25000                     | 1000pcs/bag*25bags |
| E1017   | 200–1000ul ϕ 8*71mm Fit For Finland        | 10000                     | 100pcs/bag*100bags |
| E1018   | 200–1000ul ϕ 8*71mm Fit For Switzerland    | 10000                     | 100pcs/bag*100bags |
| E1019   | 50–200ul ϕ 6*50mmFit For Switzerland       | 50000                     | 500pcs/bag*100bags |
| E1010–1 | 10ul With Filter                           | 60000                     | 1000pcs/bag*60bags |
| E1021–1 | 200ul With Filter                          | 30000                     | 1000pcs/bag*30bags |
| E1020–1 | 1000ul With Filter                         | 10000                     | 500pcs/bag*20bags  |
| E1027   | Fit for Uranus AE120                       | 9600                      | 500pcs/box*100boxs |



PIPETTE TIPS & BOXES



| Item    | Description of Goods                             | Quantity <sub>(pcs)</sub> | Packing                       |
|---------|--|---------------------------|-------------------------------|
| E1022   | Tips Box 100wells For 1000ul                     | 30                        | bulk                          |
| E1023   | Tips Box 200wells For 200ul                      | 40                        | bulk                          |
| E1024   | Tips Box 96wells For 200ul                       | 68                        | bulk                          |
| E1025   | Tips Box 96wells For 200ul&10ul                  | 68                        | bulk                          |
| E1025–1 |  |                           |                               |
| E1032   | Fit For ADDCON                                   | 10000                     | 100pcs/bag*100bags            |
| E1033   | Fit For Siemens                                  | 8640                      | 24box/carton                  |
| E1034   | Fit For TECAN Freedom EVO 200 MAC 96 Workstation | 24000                     | ( 1000pcs*2bags ) /box*12boxs |
| E1035   |  |                           |                               |
| E3100   | 1000ul Fit For TECAN                             | 4608                      | 192pcs/box*24boxs             |
| E3101   | 200ul Fit For TECAN                              | 4608                      | 192pcs/box*24boxs             |
| E4503   | 50–300ul Fit For Finland                         | 11520                     | 384pcs/box*30boxs             |




# Culture plate & pipettes


## SYSTEM

Raycaremed culture plate adopts the most advanced surface treatment technology and manufacturing process of low-temperature polymer material at present in the world. A variety of series product such as cell culture flasks, cell culture plates, isotope flasks and coplin jars can be produced to meet the needs of single cell to large-scale cell and tissue cultures in a 100 thousand stage purification plant.

## CULTURE PLATE & PIPETTES



| Item    | Description of Goods       | Quantity <sub>(pcs)</sub> | Packing                                 |
|---------|----------------------------|---------------------------|---|
| F1001   | 72 wells                   | 800                       | 1pc*800bags                             |
| F1002   | 40 wells                   | 550/600                   | 1pc*600bags with cap/1pc*550bags no cap |
| F1003-U | 96wells plate("U"bottom)   | 270                       | 1pc*270bags                             |
| F1003-V | 96wells plate("V"bottom)   | 270                       | 1pc*270bags                             |
| F1003-F | 96wells plate(Flat bottom) | 270                       | 1pc*270bags                             |
| F1004   | 24wells                    | 170                       | 1pc*170bags                             |






CULTURE PLATE & PIPETTES



| Item    | Description of Goods                                    | Quantity(pcs) | Packing     |
|---------|---|---------------|-------------|
| F1005   | 96wells detachable plat, suitable for 12 wells strip    | 270           | 1pc*270bags |
| F1006   | 52 Wells reaction plate                                 | 260           | 1pc*260bags |
| F1007   | 32 Wells rezction plate                                 | 360           | 1pc*360bags |
| F1008   | 20 Wells plate  | 180           | 1pc*180bags |
| F1021   | 96wells detachable plate, suitable for 8/12 wells strip | 270           | 1pc*270bags |
| PET250  | PET Cell Bottle 250ml                                   | 100           | 1pc*100bags |
| PET500  | PET Cell Bottle 500ml                                   | 100           | 1pc*100bags |
| PET1000 | PET Cell Bottle 1000ml                                  | 48            | 1pc*48bags  |




CULTURE PLATE & PIPETTES






Embedding box is made of durable and porous molded plastic, which can secure the samples stored in liquid and completely avoid the chemical reaction of histological solvent and provide a variety of different colors for easily identifying different tissue samples.

| Item  | Description of Goods     | Quantity(pcs) | Packing       |
|-------|--------------------------|---------------|---------------|
| F1009 | Embedding Cassette       | 5000          | 500pcs*10bags |
| F1010 | Embedding Cassette White | 4000          | 500pcs*8bags  |
| F1011 | Embedding Cassette Green | 2700          | 300pcs*9bags  |
| F1030 | Glass slides box         | 350           | 350pcs        |





TRANSFER PIPETTES

Plastic pipette fills the fragile defect of glass tube and ends the history of glass tubes and rubber plug for the transfer of high-risk contagiousness liquids. The interior of the pipette is smooth and reduces the damage to the tissue structure and cell in the liquid. It can transfer the liquid sample safely, rapidly and quantitatively. It is resistant to low temperature and can be sterilized at high temperature.



TRANSFER PIPETTES




| Item     | Description of Goods  | Quantity <sub>(pcs)</sub> | Packing             |
|----------|---|---------------------------|---------------------|
| F1014    | 25 ul & 35ul Transfer pipette                                     | 16000                     | 2000pcs*8bags       |
| F1015    | 1ml Transfer pipette  | 5000                      | 250pcs*20bags       |
| F1016    | 3ml Transfer pipette  | 5000                      | 250pcs*20bags       |
| F1016-i  | 3ml Transfer pipette paper plastic bag individual packing steaile | 2500                      | 1pcs*100bags*25bags |
| F1016-1  | 1.5ml Transfer pipette  | 5000                      | 250pcs*20bags       |
| F1016-2  | 1ml Transfer pipette  | 5000                      | 250pcs*20bags       |
| F1016-3  | 3.5ml Transfer pipette  | 4000                      | 100pcs*40bags       |
| F1016-4  | 1ml Transfer pipette  | 5000                      | 250pcs*20bags       |
| F1016-5  | 3.5ml Transfer pipette  | 4000                      | 250pcs*16bags       |
| F1016-6  | 2ml Transfer pipette  | 4000                      | 100pcs*40bags       |
| F1016-7  | 5ml Transfer pipette  | 4000                      | 250pcs*16bags       |
| F1016-8  | 1ml Transfer pipette  | 5000                      | 250pcs*20bags       |
| F1016-9  | 2ml Transfer pipette  | 5000                      | 250pcs*20bags       |
| F1016-11 | Spherical pipette   | 4000                      | 200pcs*20bags       |

# Test tubes & Racks


SYSTEM

Raycaremed test tube series is made of material PP、PS. It is transparent, smooth and high chemically tolerant. Among them, the test tube made of material PP is resistant to weak acid, weak base, organic solvent and high temperature, which can offer various specifications, flat bottom, round bottom, and various lengths.


## TEST TUBES & RACKS




| Item  | Description of Goods | Quantity <sub>(pcs)</sub> | Packing       |
|-------|----------------------|---------------------------|---------------|
| G1030 | 16*100mm             | 1600                      | 200pcs*8bags  |
| G1001 | 15*100mm, 10ml       | 2000                      | 200pcs*10bags |
| G1002 | 13*100mm, 7ml        | 3000                      | 300pcs*10bags |
| G1003 | 13*78mm, 5.5ml       | 4008                      | 334pcs*12bags |
| G1004 | 12*60mm, 3.5ml       | 6000                      | 400pcs*15bags |
| G1005 | 12*75mm, 5ml         | 4800                      | 400pcs*12bags |
| G1006 | 12*78mm, 5.5ml       | 4800                      | 400pcs*12bags |
| G1007 | 16*53mm, 6mlwith cap | 3000                      | 300pcs*10bags |




TEST TUBE & RACKS



| Item  | Description of Goods                       | Quantity <sub>(pcs)</sub> | Packing       |
|-------|--|---------------------------|---------------|
| G1009 | 13*100mm, 7ml with cap,with 5ml graduation | 3000                      | 300pcs*10bags |
| G1010 | 12*75mm, 5ml with cap,with 2ml graduation  | 4800                      | 400pcs*12bags |
| G1011 | 12*60mm, 3.5ml,with 2ml graduation         | 6000                      | 400pcs*15bags |
| G1012 | 12*75mm, 5ml,use for BD flow cytometry     | 4008                      | 334pcs*12bags |
| G1013 | 12*75mm, 5ml,hexagonal bottom              | 4800                      | 400pcs*12bags |
| G1032 | 12ml Test tube wih rubber                  | 4800                      | 400pcs*12bags |
| G1031 | 16*95mm test tube with pipette             | 1500                      | 150pcs*10bags |



TEST TUBES & RACKS




- **G1017** Equipped with the card cover. The scale is accurate to 1/4, 1/2, 2ml. The total capacity is 12ml. It can withstand 3400 RPM of centrifugal.
- The test tube mouth is provided with a aqueduct for easy drip irrigation of specimens.
- This tube is suitable for pneumatic transmission system.



| Item    | Description of Goods               | Quantity <sub>(pcs)</sub> | Packing       |
|---------|------------------------------------|---------------------------|---------------|
| G1017   | 12ml urine cup                     | 1500                      | 150pcs*10bags |
| G1017-1 | 12ml urine cup with cap            | 1500                      | 150pcs*10bags |
| G1017-2 | US urine test tue                  | 1200                      | 200pcs*6bags  |
| G1017-3 | urine test tube(new type) with cap | 2000                      | 100pcs*20bags |




TEST TUBES & RACKS




| Item  | Description of Goods  | Quantity <sub>(pcs)</sub> | Packing       |
|-------|---|---------------------------|---------------|
| G1019 | 15*50mm, 4.5ml  | 3600                      | 400pcs*9bags  |
| G1020 | 13*78mm,5.5ml   | 4800                      | 400pcs*12bags |
| G1021 | 12*75mm, 5ml PP   | 4800                      | 400pcs*12bags |
| G1027 | 15*150mm 15ml   | 1500                      | 100pcs*15bags |
| G1028 | Test tube with graduation,<br>14.7X101mm with cup.<br>Color: transparent or amber | 1000                      | 100pcs*10bags |



TEST TUBES & RACKS




| Item  | Description of Goods             | Quantity <sub>(pcs)</sub> | Packing         |
|-------|----------------------------------|---------------------------|-----------------|
| H1001 | 40Wells $\phi$ 18mm Tube Rack    | 40                        | 2pcs/bag*20bags |
| H1002 | 50Wells $\phi$ 15mm Tube Rack    | 40                        | 2pcs/bag*20bags |
| H1003 | 50Wells $\phi$ 13mm Tube Rack    | 40                        | 2pcs/bag*20bags |
| H1004 | Pillar-like $\phi$ 5*16Tube Rack | 56                        | 1pc*56bags      |
| H1005 | Pillar-like $\phi$ 5*10Tube Rack | 56                        | 1pc*56bags      |





● **H1018** 独特的底层装有硅胶垫片，采用特殊塑料做成，使试管架坚固耐用，不易破碎，能装大小各类试管。

| Item  | Description of Goods            | Quantity <sub>(pcs)</sub> | Packing    |
|-------|---------------------------------|---------------------------|------------|
| H1018 | Multi-Function Rack             | 50                        | 1pc*50bags |
| H1019 | Doble-Deek Multi -function Rack | 25                        | 1pc*25bags |

TEST TUBES & RACKS



| Item  | Description of Goods                  | Quantity (pcs) | Packing            |
|-------|---------------------------------------|----------------|--------------------|
| H1008 | Washing plate Natural White           | 200            | 1pc*200bags        |
| H1009 | 96 wells sample box for blood storage | 210            | in bulk            |
| F1019 | Reservoir type "V"                    | 350            | Individual Packing |
| E1026 | Pipette shelf                         | 10             | in bulk            |



Cuvettes

SYSTEM

Raycaremed colorimetric cup series adopts imported optical plastics, and is casted with high precision mould. The inside and outside surface of the product is transparent, level and smooth with superior light transmission performance. It is suitable for all kinds of automatic biochemical analyzer, blood coagulation instrument and spectrophotometer at home and abroad.


CUVETTES




| Item    | Description of Goods                  | Quantity (pcs) | Packing                |
|---------|---------------------------------------|----------------|------------------------|
| A1001   | Tianjin Pulisheng & French stago      | 1600           | 50pcs/box*32boxes/ctn  |
| A1001-1 | Steellex                              | 2300           | 72pcs/bag*32bags-4pcs  |
| A1002-1 | American MD 3000 Biochemical Analysis | 960            | 240pcs/box*4boxes/ctn  |
| A1003   | CA100 mid-Asia                        | 15000          | 1000pcs/bag*15bags/ctn |
| A1004   | French Organon                        | 2688           | 84pcs/bag*32bags/ctn   |
| A1005   | BN100 Special Protein Apparatus       | 2688           | 70pcs/box*32boxes/ctn  |



CUVETTES



| Item    | Description of Goods           | Quantity (pcs) | Packing                |
|---------|--------------------------------|----------------|------------------------|
| A1006   | Roche Biochemical Analysis Cup | 500            | 1pc*500bags            |
| A1007   | American Pacific, German TECO  | 3200           | 100pcs/bag*32bags/ctn  |
| A1007-1 | MC Coatron 1800                | 3200           | 100pcs/bag*32bags/ctn  |
| A1008   | Beckman                        | 1500           | 100pcs/box*15boxes/ctn |
| A1008-1 | Beckman TOP                    | 1500           | 100pcs/box*15boxes/ctn |




CUVETTES



| Item  | Description of Goods          | Quantity(pcs) | Packing                           |
|-------|-------------------------------|---------------|-----------------------------------|
| A1010 | Finland Kangyi 12 Holes       | 3200          | 20pcs*5boxes/bag*32bags           |
| A1011 | American Corning cups 5 Holes | 2366          | 13pcs/box*13boxes/bag*14boxes/ctn |
| A1012 | Germany BECL4 8Holes          | 1000          | 25pcs/box*10boxes/bag*4boxes/ctn  |
| A1015 | German BE 4 Holes             | 3500          | 500pcs*7bags                      |
| A1017 | German BE 1 Hole              | 24000         | 2000pcs*12bags                    |
| A1018 | Spain Cruor Apparatus         | 8000          | 1000pcs*8bags                     |
| A1019 | German LABOR                  | 10000         | 1000pcs*10bags                    |
| A1020 | German TECO Sample Cup        | 12000         | 1000pcs*12bags                    |



CUVETTES



| Item  | Description of Goods     | Quantity(pcs) | Packing         |
|-------|--------------------------|---------------|-----------------|
| B1001 | German BE Automatic      | 2610          | 261pcs*10bags   |
| B1002 | German SIGMA KC-40       | 12000         | 1000pcs*12bags  |
| B1005 | Yapei AXSYM              | 3600          | 400pcs*9bags    |
| B1006 | Yapei AXSYM reaction jar | 1200          | 200pcs*6bags    |
| B1007 | Blood tube               | 3000          | 300pcs*10bags   |
| B1009 | Japan Dongya CA1500      | 20000         | 2000pcs*10bags  |
| B1010 | Japan Dongya CA530       | 20000         | 2000pcs*10bags  |
| B1011 | Japan Dongya Ca6000      | 20000         | 2000pcs*10bags  |
| B1020 | CA7000                   | 20000         | 1000pcs*20boxes |





CUVETTES





| Item  | Description of Goods                 | Quantity (pcs) | Packing         |
|-------|--------------------------------------|----------------|-----------------|
| B1025 | Yapei R2000                          | 10000          | 1000pcs*10boxes |
| B1026 | Siemens ADVIA Centaur CP             | 6000           | 1000pcs*6bags   |
| B1027 | Siphon quantitative dropping pipette | 10000          | 1000pcs*10bags  |
| B1015 | Rayto RAC-1800 hematology machine    | 20000          | 1000pcs*20boxes |



CUVETTES




| Item  | Description of Goods                                   | Quantity (pcs) | Packing           |
|-------|--|----------------|-------------------|
| B1012 | American Beckman SpecialProtein<br>ARRY-360 Sample Cup | 400            | 40pcs*10bags      |
| B1013 | Color Dilution Plate                                   | 100            | 1*100box          |
| B1014 | Baolingman 10holes                                     | 1000           | 200pcs*5boxes     |
| B1016 | For Special Protein Analyzer                           | 960            | 40pcs/box*24boxes |
| B1018 | Laibang  | 20000          | 2000pcs*10bags    |



CUVETTES



| Item  | Description of Goods  | Quantity <sub>(pcs)</sub> | Packing        |
|-------|---|---------------------------|----------------|
| C1001 | Hitachi 7150,7060 Series;Aulinbas 2ml                                       | 5000                      | 500pcs*10bags  |
| C1002 | American Beckman CX Series,0.5ml  | 8000                      | 1000pcs*8bags  |
| C1003 | IL Sample Cup,Monarch Biochemical Analysis Apparatus(0.5ml)                 | 13000                     | 1000pcs*13bags |
| C1004 | American Beckman 700(1.5ml), American Dubang, Holland<br>VITALAB SAMPLE CUP | 8000                      | 1000pcs*8bags  |
| C1025 | FC serum sample cup (L)   | 10000                     | 1000pcs*10bags |
| C1026 | FC serum sample cup (S)   | 10000                     | 1000pcs*10bags |
| C1028 | Beckman DXI800RV  | 12000                     | 1000pcs*12bags |
| C1029 | Analysis tube   | 2400                      | 200pcs*12bags  |




CUVETTES



| Item  | Description of Goods           | Quantity <sub>(pcs)</sub> | Packing         |
|-------|--------------------------------|---------------------------|-----------------|
| C1006 | AC900(22*45)                   | 2100                      | 100pcs*21bags   |
| C1010 | Mairui(30*50)                  | 10000                     | 65pcs*153.8bags |
| C1016 | F800 100holes sample rack      | 11                        | 1pc*11bags      |
| C1017 | Hitachi Minimum                | 6000                      | 500pcs*12bags   |
| C1019 | Hitachi Conical                | 9000                      | 1000pcs*9bags   |
| C1023 | An'Tai 30 Holes Reaction Plate | 900                       | 1pc*900bags     |
| C1024 | Antai sample cup               | 2600                      | 2000pcs*13bags  |



CUVETTES



| Item    | Description of Goods                            | Quantity(pcs) | Packing          |
|---------|---|---------------|------------------|
| C1011   | Holl ACE  | 24000         | 2000pcs*12bags   |
| C1013   | American M10,Italy BT2245 Biochemical Apparatus | 20000         | 2000pcs*10bags   |
| C1014   | American Kangning 560, 550 Series               | 18000         | 2000pcs*9bags    |
| C1015   | American M10,Italy BT2245 Biochemical Apparatus | 10000         | 1000pcs*10bags   |
| C1020   | Kodak Sample Cup                                | 12000         | 1000pcs*12bags   |
| C1021   | Baolingman Sample Tube                          | 40000         | 5000pcs*8bags    |
| C1022   | Kodak tip                                       | 16000         | 1000pcs*40bags   |
| C1021-1 | Baolingman black pipette tips                   | 20000         | 1000pcs*20boxes  |
| C1021-2 | Roche Tip&Cup                                   | 48bo          | 6boxes/bag*8bags |
| C1021-5 | Roche Tip(with box)                             | 120           | 120boxes         |
| C1021-6 | Roche Cup(with box)                             | 240           | 240boxes         |



PHARMACEUTICAL PACKAGING MATERIAL

